

## Product datasheet for **RC201834L2V**

### **HNRPH1 (HNRNPH1) (NM\_005520) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HNRPH1 (HNRNPH1) (NM_005520) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HNRPH1
Synonyms:	hnRNPH; HNRPH; HNRPH1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_005520
ORF Size:	1347 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201834).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005520.1</a>
RefSeq Size:	2274 bp
RefSeq ORF:	1350 bp
Locus ID:	3187
UniProt ID:	<a href="#">P31943</a>
Cytogenetics:	5q35.3
Domains:	RRM
MW:	49.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of a subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA. These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some may shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNA and is very similar to the family member HNRPF. This gene may be associated with hereditary lymphedema type I. Alternatively spliced transcript variants have been described [provided by RefSeq, Mar 2012]