

Product datasheet for **RC201861A1V**

Human RNF141 (NM_016422) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RNF141 (NM_016422) AAV Particle
Tag: Myc-DDK
Symbol: RNF141
Synonyms: RFP141; ZFP26; ZNF230
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201861 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGACAGCAAATTTTCGGATCAGACACAGTTGGTTATTAACAAGTTACCAGAAAAAGTAGCAAAACATG
 TTACGTTGGTTCGAGAGAGTGGCTCCTTAACCTATGAAGAATTTCTCGGGAGAGTAGCTGAGCTTAATGA
 TGTAACGGCTAAAGTGGCTTCTGGCCAGGAAAAACATCTTCTCTTTGAGGTACAACCTGGGTCTGATTCC
 TCTGCTTTTGGAAAGTGGTGTACGGGTGGTCTGTACCAAGATTAACAAAAGCAGTGGCATTGTGGAGG
 CATCACGGATCATGAATTTATACAGTTTATTCAACTTTATAAAGATATCACAAGTCAAGCAGCAGGAGT
 ATTGGCACAGAGCTCCACCTCTGAAGAACCTGATGAAAACCTCATCCTCTGTAAACATCTTGTGAGGCTAGT
 CTTTGGATGGGAAGGGTGAAGCAGCTGACCGATGAGGAGGAGTGTGTATCTGTATGGATGGGCGGGCTG
 ACCTCATCCTGCCTTGTGCTCACAGCTTTTGTGAGAAGTGTATTGATAAATGGAGTGATCGACACAGGAA
 TTGCCCTATTTGTGCGCTACAGATGACTGGAGCAAAATGAATCTTGGGTGGTATCAGATGCACCCACTGAA
 GATGATATGGCTAACTATATTCTTAACATGGCTGATGAGGCAGGCCAGCCCCACAGGCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201861 protein sequence
 Red=Cloning site Green=Tags(s)

MGQQISDQTQLVINKLPEKVAKHVTLVRESGLTYEEFLGRVAELNDVTAKVASGQEKHLLFEVQPGSDS
 SAFWKVVVRVCTKINKSSGIVEASRIMNLYQFIQLYKDITSQAAGVLAQSSTSEEPDENSSSVTSCQAS
 LWMGRVKQLTDEEECCICMDGRADLILPCAHSFCQKCIDKWSDRHRNCPICRLQMTGANESWVVSADPTE
 DDMANYILNMADEAGQPHRP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_016422
ORF Size:	690 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_016422.3
RefSeq Size:	4084 bp
RefSeq ORF:	693 bp
Locus ID:	50862
UniProt ID:	Q8WVD5
Cytogenetics:	11p15
MW:	25.5 kDa