

## Product datasheet for **RC201868L4V**

### CGI 62 (ZC2HC1A) (NM\_016010) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CGI 62 (ZC2HC1A) (NM_016010) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CGI 62
Synonyms:	C8orf70; CGI-62; FAM164A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016010
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201868).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016010.1</a>
RefSeq Size:	3352 bp
RefSeq ORF:	978 bp
Locus ID:	51101
UniProt ID:	<a href="#">Q96GY0</a>
Cytogenetics:	8q21.13
Domains:	zf-C2H2
MW:	35.1 kDa


[View online »](#)