

## Product datasheet for **RC201869L2V**

### ABHD5 (NM\_016006) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ABHD5 (NM_016006) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABHD5
Synonyms:	CGI58; IECN2; NCIE2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_016006
ORF Size:	1047 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201869).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016006.3</a>
RefSeq Size:	5370 bp
RefSeq ORF:	1050 bp
Locus ID:	51099
UniProt ID:	<a href="#">Q8WTS1</a>
Cytogenetics:	3p21.33
Domains:	abhydrolase
Protein Families:	Protease



[View online »](#)

**MW:** 39.1 kDa

**Gene Summary:** The protein encoded by this gene belongs to a large family of proteins defined by an alpha/beta hydrolase fold, and contains three sequence motifs that correspond to a catalytic triad found in the esterase/lipase/thioesterase subfamily. It differs from other members of this subfamily in that its putative catalytic triad contains an asparagine instead of the serine residue. Mutations in this gene have been associated with Chanarin-Dorfman syndrome, a triglyceride storage disease with impaired long-chain fatty acid oxidation. [provided by RefSeq, Jul 2008]