

## Product datasheet for **RC201888L3V**

### TDP1 (NM\_001008744) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TDP1 (NM_001008744) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TDP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001008744
ORF Size:	1824 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201888).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001008744.1</a>
RefSeq Size:	3540 bp
RefSeq ORF:	1827 bp
Locus ID:	55775
UniProt ID:	<a href="#">Q9NUW8</a>
Cytogenetics:	14q32.11
Protein Families:	Druggable Genome
MW:	68.4 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is involved in repairing stalled topoisomerase I-DNA complexes by catalyzing the hydrolysis of the phosphodiester bond between the tyrosine residue of topoisomerase I and the 3-prime phosphate of DNA. This protein may also remove glycolate from single-stranded DNA containing 3-prime phosphoglycolate, suggesting a role in repair of free-radical mediated DNA double-strand breaks. This gene is a member of the phospholipase D family and contains two PLD phosphodiesterase domains. Mutations in this gene are associated with the disease spinocerebellar ataxia with axonal neuropathy (SCAN1). [provided by RefSeq, Aug 2016]