

Product datasheet for **RC201893A1V**

Human CLIC4 (NM_013943) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CLIC4 (NM_013943) AAV Particle
Tag:	Myc-DDK
Symbol:	CLIC4
Synonyms:	CLIC4L; H1; huH1; MTCLIC; p64H1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201893 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGTTGTCGATGCCGCTGAATGGGCTGAAGGAGGAGGACAAAGAGCCCTCATCGAGCTCTTCGTCA
AGGCTGGCAGTGATGGTGAAAGCATAGGAACTGCCCTTTTCCAGAGGCTCTTCATGATTCTTTGGCT
CAAAGGAGTTGTATTTAGTGTGACGACTGTTGACCTGAAAAGGAAGCCAGCAGACCTGCAGAACTTGGCT
CCCGGGACCCACCCACCATTTATACTTTCAACAGTGAAGTCAAAACGGATGTAATAAGATTGAGGAAT
TTCTTGAAGAAGTCTTATGCCCTCCCAAGTACTTAAAGCTTTCACCAAAACACCCAGAATCAAATACTGC
TGGAATGGACATCTTTGCCAAATTCCTGTCATATATCAAGAATTCAAGGCCAGAGGCTAATGAAGCACTG
GAGAGGGGTCTCCTGAAAACCTGCAGAACTGGATGAATATCTGAATTCTCCTCTCCCTGATGAAATTG
ATGAAAATAGTATGGAGGACATAAAGTTTTCTACACGTAATTTCTGGATGGCAATGAAATGACATTAGC
TGATTGCAACCTGCTGCCAACTGCATATTGTCAAGGTGGTGCCAAAAAATATCGCAACTTTGATATT
CCAAAAGAAATGACTGGCATCTGGAGATACCTAATAATGCATACAGTAGGGACGAGTTCACCAATACCT
GTCCAGTGATAAGGAGGTTGAAATAGCATATAGTGATGTAGCCAAAAGACTCACCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201893 protein sequence Red=Cloning site Green=Tags(s) MALSMPLNGLKEEDKEPLIELFVKAGSDGESIGNCPFSQRLFMILWLKGVVFSVTTVDLKRKPADLQNLAPGTHPPFITFNSEVKTDVNKIEEFLEEVLCPPKYLKLSPKHPESNTAGMDIFAKFSAYIKNSRPEANEALERGLLKTQLDEYLNPLPDEIDENSMEDIKFSTRKFLDGNEMTLADCNLLPKLHIVKVVAKKYRNFDIPKEMTGIWRYLTNAYSRDEFTNTCPSDKEVEIAYSDEVAKRLTK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_013943
ORF Size:	759 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_013943.1
RefSeq Size:	4452 bp
RefSeq ORF:	762 bp
Locus ID:	25932
UniProt ID:	Q9Y696
Cytogenetics:	1p36.11
MW:	28.8 kDa