

Product datasheet for **RC201894A1V**

Human HMG20A (NM_018200) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HMG20A (NM_018200) AAV Particle
Tag:	Myc-DDK
Symbol:	HMG20A
Synonyms:	HMGX1; HMGXB1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201894 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAAACCTTGATGACTAGCTCCACCCTACCGCCCTTTTTCAGATGAAGACGGTTCCAAGGAGAGTA
ATGATCTGGCTACCACTGGGTAAATCACCCAGAGGTTCCATACAGTAGTGGCGCCACATCATCCACCAA
CAATCCAGAATTTGTGGAGGATCTCTCAAGGTCAGTTGCTTCAGAGTGAGTCTTCAAATGCAGCAGAA
GGCAATGAACAGAGGCATGAAGATGAGCAACGAAGTAACGAGGAGGTTGGTCCAAAGGAAGAAAGAGGA
AGAAACCTCTTCGAGACAGCAATGCACCCAAATCCCCCTTACAGGATATGTTTCGGTTCATGAATGAGCG
TCGAGAAACAACCTTCGAGCAAAAGAGACCAGAAAGTCCCATTTCCAGAAATCACAAGGATGTTAGGCAATGAA
TGGAGTAACTGCCTCCTGAGGAAAAACAGCGCTACCTTGATGAAGCAGACAGAGATAAGGAGCGTTACA
TGAAGGAACCTGGAACAGTATCAGAAAACAGAGGCTACAAGGTCTTCAGTAGGAAAACCCAGGACCGTCA
GAAAGGCAATCTCATAGGCAAGATGCAGCCCGCAGGCCACTCATGATCATGAGAAAGAAACAGAGGTA
AAGGAACGGTCTGTTTTTGACATCCCTATATTTACAGAGGAATCTTGAACCATAGCAAAGCTCGGGAAG
CAGAGCTCCGCCAGCTTCGAAATCCAACATGGAGTTTGAGGAGAGGAATGCAGCCCTGCAAAAGCACGT
GGAGAGCATGCGCACAGCAGTGGAGAAGCTGGAGGTGGATGTATCCAGGAGCGGAGCCGCAACACAGTC
TTACAGCAGCACCTGGAGACCCTGCGGCAGGTGCTGACCAGCAGCTTTGCCAGCATGCCCTTGCTGGAA
GTGGAGAGACACCTACAGTGGACACCATTTGACTCATATATGAACAGACTGCACAGTATTATTTTAGCTAA
TCCCAAGACAATGAAAACCTTACATAGCTACAGTTTCGAGAAGTTGTGAACAGACTCGATCGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201894 protein sequence Red=Cloning site Green=Tags(s) MENLMTSSTLPPLFADEDGSKESNDLATTGLNHPEVPYSSGATSSSTNNPEFVEDLSQGQLLQSESSNAAE GNEQRHEDEQRSKRGGWSKGRKRKKPLRDSNAPKSPLTGYVRFMNERREQLRAKRPEVPFPEITRMLGNE WSKLPPEEKQRYLDEADRDKERYMKELEQYQKTEAYKVF SRKTQDRQKGKSHRQDAARQATHDHEKETEV KERSVFDIPIFTEEFNLHNSKAREELRQLRKS NMEFEERNAALQKHVESMRTAVEKLEVDVIQERSRNTV LQQHLET LRQVLTSSFASMP LPGSGETPTVDTIDS YMNRLHSIILANPQDNENFIATVREV VNR LDR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018200
ORF Size:	1041 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_018200.2</u>
RefSeq Size:	3898 bp
RefSeq ORF:	1044 bp
Locus ID:	10363
UniProt ID:	<u>Q9NP66</u>
Cytogenetics:	15q24.3
MW:	40.1 kDa