

Product datasheet for **RC201894L4V**

HMG20A (NM_018200) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HMG20A (NM_018200) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HMG20A
Synonyms:	HMGX1; HMGXB1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018200
ORF Size:	1041 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201894).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018200.2
RefSeq Size:	3898 bp
RefSeq ORF:	1044 bp
Locus ID:	10363
UniProt ID:	Q9NP66
Cytogenetics:	15q24.3
Domains:	HMG
Protein Families:	Transcription Factors



[View online »](#)

MW: 40.1 kDa

Gene Summary: Plays a role in neuronal differentiation as chromatin-associated protein. Acts as inhibitor of HMG20B. Overcomes the repressive effects of the neuronal silencer REST and induces the activation of neuronal-specific genes. Involved in the recruitment of the histone methyltransferase KMT2A/MLL1 and consequent increased methylation of histone H3 lysine 4 (By similarity).[UniProtKB/Swiss-Prot Function]