

Product datasheet for **RC201901A1V**

Human ASB9 (NM_001031739) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human ASB9 (NM_001031739) AAV Particle

Tag: Myc-DDK

Symbol: ASB9

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide

Sequence: >RC201901 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGATGGCAAACAAGGGGGCATGGATGGGAGCAAGCCCGGGGGCCAAGGGACTTTCTGGCATCAGGC
TTCTTTCAAACCCATTGATGGGCGATGCTGTGTCTGATTGGTCTCCTATGCATGAAGCTGCAATCCACGG
ACATCAGCTGTCTCTGAGGAACCTCATCAGCCAGGGGTGGGCTGTGAACATCATCACGGCAGATCATGTT
TCCCCACTCCATGAAGCCTGTCTTGGAGGTCATCTCTTGTGTGAAGATTTTATTAAGCATGGAGCTC
AGGTGAATGGTGTGACAGCAGACTGGCACACTCCACTGTTTAATGCTTGTGTGACGGCAGCTGGGATTG
TGTGAATTTGCTTCTGCAGCACGGAGCCAGCGTTCAACCTGAGAGTGATCTGGCATCCCCATCCATGAA
GCTGCTAGGAGAGGCCACGTGGAGTGTGTCAACTCTCTTATAGCTTATGGGGGCAACATTGACCATAAGA
TCAGCCACCTGGGCACTCCACTCTATTTGGCTTGTGAAAACCAACAGAGAGCCTGTGTCAAGAAGCTTCT
GGAGTCAGGAGCGGACGTGAACCAAGGAAAGGTGAGGATCCCCACTTCATGCAGTGGCCAGGACAGCC
AGTGAAGAGCTGGCCTGCCTGCTCATGGATTTTGGAGCGGACACCCAGGCCAAGAATGCTGAAGGCAAAC
GTCCTGTGGAGCTGGTGCCTCCAGAGAGCCCTTGGCCCAGCTCTTCTTGGAGAGAGAAGGGCCCCCTTC
TTTGATGCAGTTATGCCGCTTAGAATTCGGAAGTGTGTTGGAATCCAGCAGCATCATAAGATAACCAAA
CTCGTCCTCCAGAGGATCTGAAACAGTTTCTCTACATCTT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201901 protein sequence Red=Cloning site Green=Tags(s) MDGKQGGMDGSKPAGPRDFPGIRLLSNPLMGDAVSDWSPMHEAAIHGHQLSLRNLIISQGWAVNIITADHV SPLHEACLGGLSCVKILLKHGAQVNGVTADWHTPLFNACVSGSWDCVNLLQHGASVQPESDLASPIHE AARRGHVECVNSLIAYGGNIDHKISHLGTPLYLACENQQRACVKKLLESGADVNQGGQDSPLHAVARTA SEELACLLMDFGADTQAKNAEGKRPVELVPPEPLAQLFLEREGPPSLMQLCRLRIRKCFGIQHHKITS LVLPEDLKQFLLHL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001031739
ORF Size:	882 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001031739.1</u>
RefSeq Size:	1714 bp
RefSeq ORF:	885 bp
Locus ID:	140462
UniProt ID:	<u>Q96DX5</u>
Cytogenetics:	Xp22.2
MW:	31.9 kDa