

Product datasheet for **RC201919L4V**

HDAC11 (NM_024827) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HDAC11 (NM_024827) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HDAC11
Synonyms:	HD11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024827
ORF Size:	1041 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201919).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024827.1
RefSeq Size:	2918 bp
RefSeq ORF:	1044 bp
Locus ID:	79885
UniProt ID:	Q96DB2
Cytogenetics:	3p25.1
Protein Families:	Druggable Genome, Transcription Factors
MW:	39.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a class IV histone deacetylase. The encoded protein is localized to the nucleus and may be involved in regulating the expression of interleukin 10. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Apr 2009]