

Product datasheet for **RC201924A1V**

Human Calretinin (CALB2) (NM_001740) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Calretinin (CALB2) (NM_001740) AAV Particle
Tag:	Myc-DDK
Symbol:	Calretinin
Synonyms:	CAB29; CAL2; CR
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201924 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGCCCGCAGCAGCAGCCCCCTTACCTGCACCTGGCCGAGCTGACGGCGTCCCAGTTCCTGGAAA
TATGGAAGCACTTTGACGCAGACGAAATGGGTATATTGAAGGTAAAGAGCTAGAAAACCTTTTCCAAGA
GCTGGAGAAGGCAAGGAAAGGCTCTGGCATGATGTCAAAGAGTGACAACCTTTGGAGAAAAGATGAAGGAG
TTCATGCAGAAGTATGATAAACTCAGATGGGAAATCGAGATGGCAGAGCTGGCGCAGATCCTGCCAA
CCGAAGAGAATTCTTCTGTGCTTCAGGCAGCACGTGGGCTCCAGCACCGAGTTTATGGAGGCTTGGCG
GAAGTACGACACAGACAGGAGTGGCTACATCGAAGCCAATGAGCTCAAGGGATTCTGTGACACCTGCTG
AAGAAGGCGAACC GGCGTACGATGAGCCCAAGCTCCAGGAATACACCCAAACCATACTACGGATGTTTG
ACTTGAACGGGGATGGCAAATTGGGCTCTCAGAGATGTCCCGACTCCTGCCTGTCCAGGAAAACCTCCT
GCTTAAATTTAGGGCATGAAGCTGACCTCAGAGGAGTTTAACGCGATCTTCACATTTTACGACAAGGAT
AGAAGCGGCTACATTGACGAGCATGAGCTGGATGCCCTTTTGAAGGATCTGTACGAGAAAAACAAAAGG
AAATGAATATTCAACAGCTACCAACTACAGAAAGAGCGTCATGTCCTTGGCAGAGGCAGGAAGCTCTA
CCGCAAGGACCTGGAGATTGTGCTCTGCAGCGAGCCCCCATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201924 protein sequence Red=Cloning site Green=Tags(s) MAGPQQQPPYLHLAELTASQFLEIWKHFDADGNGYIEGKELENFFQELEKARKGSGMMSKSDNFGEKMKE FMQKYDKNSDGKIEMAELAQILPTEENFLLCFRQHVGSSTEFMEAWRKYDTRSGYIEANELKGFLSDLL KKANRPYDEPKLQEYTQTILRMFDLNGDGKLGLEMSRLLPVQENFLLKFQGMKLTSEEFNAIFTFYDKD RSGYIDEHELDALLKDLYEKNKKEMNIQQLTNYRKSVMSLAEAGKLYRKDLEIVLCSEPPM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001740
ORF Size:	813 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001740.2</u>
RefSeq Size:	1485 bp
RefSeq ORF:	816 bp
Locus ID:	794
UniProt ID:	<u>P22676</u>
Cytogenetics:	16q22.2
MW:	31.6 kDa