

## Product datasheet for **RC201927L4V**

### TBCE (NM\_003193) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TBCE (NM_003193) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TBCE
Synonyms:	HRD; KCS; KCS1; pac2; PEAMO
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003193
ORF Size:	1581 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201927).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003193.3</a>
RefSeq Size:	1977 bp
RefSeq ORF:	1584 bp
Locus ID:	6905
UniProt ID:	<a href="#">Q15813</a>
Cytogenetics:	1q42.3
Domains:	CAP_GLY, LRR
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 59.3 kDa

**Gene Summary:** Cofactor E is one of four proteins (cofactors A, D, E, and C) involved in the pathway leading to correctly folded beta-tubulin from folding intermediates. Cofactors A and D are believed to play a role in capturing and stabilizing beta-tubulin intermediates in a quasi-native confirmation. Cofactor E binds to the cofactor D/beta-tubulin complex; interaction with cofactor C then causes the release of beta-tubulin polypeptides that are committed to the native state. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]