

## Product datasheet for **RC201930L2V**

### TRIM31 (NM\_007028) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRIM31 (NM_007028) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRIM31
Synonyms:	C6orf13; HCG1; HCGI; RNF
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_007028
ORF Size:	1275 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201930).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007028.3</a>
RefSeq Size:	2049 bp
RefSeq ORF:	1278 bp
Locus ID:	11074
UniProt ID:	<a href="#">Q9BZY9</a>
Cytogenetics:	6p22.1
MW:	48.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that functions as an E3 ubiquitin-protein ligase. This gene shows altered expression in certain tumors and may be a negative regulator of cell growth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]