

Product datasheet for **RC201931L4V**

DNAJC12 (NM_021800) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNAJC12 (NM_021800) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNAJC12
Synonyms:	HPANBH4; JDP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021800
ORF Size:	594 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201931).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021800.2
RefSeq Size:	1218 bp
RefSeq ORF:	597 bp
Locus ID:	56521
UniProt ID:	Q9UKB3
Cytogenetics:	10q21.3
Domains:	Dnaj
MW:	23.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]