

Product datasheet for **RC201935L4V**

IFI44L (NM_006820) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IFI44L (NM_006820) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IFI44L
Synonyms:	C1orf29; GS3686
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006820
ORF Size:	1239 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201935).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006820.2 , NP_006811.1
RefSeq Size:	5889 bp
RefSeq ORF:	1359 bp
Locus ID:	10964
UniProt ID:	Q53G44
Cytogenetics:	1p31.1
MW:	46.9 kDa
Gene Summary:	Exhibits a low antiviral activity against hepatitis C virus.[UniProtKB/Swiss-Prot Function]



[View online »](#)