

Product datasheet for **RC201936L4V**

PNRC1 (NM_006813) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PNRC1 (NM_006813) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PNRC1
Synonyms:	B4-2; PNAS-145; PROL2; PRR2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006813
ORF Size:	981 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201936).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006813.1
RefSeq Size:	2133 bp
RefSeq ORF:	984 bp
Locus ID:	10957
UniProt ID:	Q12796
Cytogenetics:	6q15
Protein Families:	Druggable Genome, Transcription Factors
MW:	35.2 kDa



[View online »](#)

Gene Summary: Nuclear receptor coactivator. May play a role in signal transduction.[UniProtKB/Swiss-Prot Function]