

Product datasheet for **RC201940L1V**

AHCYL1 (NM_006621) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AHCYL1 (NM_006621) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AHCYL1
Synonyms:	DCAL; IRBIT; PPP1R78; PRO0233; XPVKONA
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_006621
ORF Size:	1590 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201940).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006621.4
RefSeq Size:	3949 bp
RefSeq ORF:	1593 bp
Locus ID:	10768
UniProt ID:	O43865
Cytogenetics:	1p13.3
Domains:	AdoHcyase
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Cysteine and methionine metabolism, Metabolic pathways, Selenoamino acid metabolism

MW: 59 kDa

Gene Summary: The protein encoded by this gene interacts with inositol 1,4,5-trisphosphate receptor, type 1 and may be involved in the conversion of S-adenosyl-L-homocysteine to L-homocysteine and adenosine. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jun 2011]