

## Product datasheet for **RC201944L4V**

### TRIM38 (NM\_006355) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRIM38 (NM_006355) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRIM38
Synonyms:	RNF15; RORET
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006355
ORF Size:	1395 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201944).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006355.2</a>
RefSeq Size:	5629 bp
RefSeq ORF:	1398 bp
Locus ID:	10475
UniProt ID:	<a href="#">O00635</a>
Cytogenetics:	6p22.2
Domains:	zf-B_box, RING, SPRY, PRY
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 53.4 kDa

**Gene Summary:** This gene encodes a member of the tripartite motif (TRIM) family. The encoded protein contains a RING-type zinc finger, B box-type zinc finger and SPRY domain. The function of this protein has not been identified. A pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Jul 2012]