

## Product datasheet for **RC201946L2V**

### **SSSCA1 (ZNRD2) (NM\_006396) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SSSCA1 (ZNRD2) (NM_006396) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNRD2
Synonyms:	p27; SSSCA1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006396
ORF Size:	597 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201946).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006396.1</a>
RefSeq Size:	661 bp
RefSeq ORF:	600 bp
Locus ID:	10534
UniProt ID:	<a href="#">O60232</a>
Cytogenetics:	11q13.1
MW:	21.5 kDa



[View online »](#)

**Gene Summary:**

This antigen is recognized by a subset of anti-centromere antibodies from patients with scleroderma and/or Sjogren's syndrome. Subcellular localization has not yet been established. [provided by RefSeq, Jul 2008]