

## Product datasheet for **RC201964L4V**

### DPEP1 (NM\_004413) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DPEP1 (NM_004413) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DPEP1
Synonyms:	MBD1; MDP; RDP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004413
ORF Size:	1233 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201964).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004413.1</a>
RefSeq Size:	1752 bp
RefSeq ORF:	1236 bp
Locus ID:	1800
UniProt ID:	<a href="#">P16444</a>
Cytogenetics:	16q24.3
Domains:	Renal_dipeptase
Protein Families:	Protease



[View online »](#)

**MW:** 45.7 kDa

**Gene Summary:** The protein encoded by this gene is a kidney membrane enzyme involved in the metabolism of glutathione and other similar proteins by dipeptide hydrolysis. The encoded protein is known to regulate leukotriene activity by catalyzing the conversion of leukotriene D4 to leukotriene E4. This protein uses zinc as a cofactor and acts as a disulfide-linked homodimer. [provided by RefSeq, Dec 2020]