

Product datasheet for **RC201965A1V**

Human **BCL2A1** (NM_004049) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human BCL2A1 (NM_004049) AAV Particle
Tag: Myc-DDK
Symbol: BCL2A1
Synonyms: ACC-1; ACC-2; ACC1; ACC2; BCL2L5; BFL1; GRS; HBPA1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201965 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACAGACTGTGAATTTGGATATATTTACAGGCTGGCTCAGGACTATCTGCAGTGCGTCTACAGATAC
CACAACCTGGATCAGGTCCAAGCAAAACGTCCAGAGTGCTACAAAATGTTGCGTTCTCAGTCCAAAAAGA
AGTGGAAAAGAATCTGAAGTCATGCTTGGACAATGTTAATGTTGTGTCGGTAGACACTGCCAGAACACTA
TTCAACCAAGTGATGGAAAAGGAGTTTGAAGACGGCATCATTAACTGGGAAGAATTGTAACCATATTTG
CATTTGAAGGTATTCTCATCAAGAACTTCTACGACAGCAAATTGCCCGGATGTGGATACCTATAAGGA
GATTTTCATATTTTGTGCGGAGTTCATAATGAATAACACAGGAGAATGGATAAGGCAAAACGGAGGCTGG
GAAAATGGCTTTGTAAAGAAGTTTGAACCTAAATCTGGCTGGATGACTTTTCTAGAAGTTACAGGAAAGA
TCTGTGAAATGCTATCTCTCTGAAGCAATACTGT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201965 protein sequence
Red=Cloning site Green=Tags(s)

MTDCEFGYIYRLAQDYLCVLQIPQPGSGPSKTSRVLQNVAFSVQKEVEKNLKSCLDNVNVVSDTARTL
FNQVMEKEFEDGIINWGRIVTIFAFEGILIKLLRQQIAPDVTYKEISYFVAEFIMNNTGEWIRQNGGW
ENGFVKKFEPKSGWMTFLEVTGKICEMLSLLKQYC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_004049
ORF Size:	525 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004049.2
RefSeq Size:	899 bp
RefSeq ORF:	528 bp
Locus ID:	597
UniProt ID:	Q16548
Cytogenetics:	15q25.1
MW:	20.1 kDa