

Product datasheet for **RC201972L4V**

Methylmalonyl Coenzyme A mutase (MUT) (NM_000255) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Methylmalonyl Coenzyme A mutase (MUT) (NM_000255) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Methylmalonyl Coenzyme A mutase
Synonyms:	MCM; MUT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000255
ORF Size:	2250 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201972).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000255.1
RefSeq Size:	3886 bp
RefSeq ORF:	2253 bp
Locus ID:	4594
UniProt ID:	P22033
Cytogenetics:	6p12.3
Domains:	MM_CoA_mutase, B12-binding



[View online »](#)

Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Propanoate metabolism, Valine, leucine and isoleucine degradation
MW:	83.2 kDa
Gene Summary:	This gene encodes the mitochondrial enzyme methylmalonyl Coenzyme A mutase. In humans, the product of this gene is a vitamin B12-dependent enzyme which catalyzes the isomerization of methylmalonyl-CoA to succinyl-CoA, while in other species this enzyme may have different functions. Mutations in this gene may lead to various types of methylmalonic aciduria. [provided by RefSeq, Jul 2008]