

Product datasheet for **RC201973L1V**

Fatty Acid Binding Protein 5 (FABP5) (NM_001444) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fatty Acid Binding Protein 5 (FABP5) (NM_001444) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Fatty Acid Binding Protein 5
Synonyms:	E-FABP; EFABP; KFABP; PA-FABP; PAFABP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001444
ORF Size:	405 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201973).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001444.1
RefSeq Size:	751 bp
RefSeq ORF:	408 bp
Locus ID:	2171
UniProt ID:	Q01469
Cytogenetics:	8q21.13
Domains:	lipocalin



[View online »](#)

Protein Pathways: PPAR signaling pathway

MW: 15.2 kDa

Gene Summary: This gene encodes the fatty acid binding protein found in epidermal cells, and was first identified as being upregulated in psoriasis tissue. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. FABPs may play roles in fatty acid uptake, transport, and metabolism. Polymorphisms in this gene are associated with type 2 diabetes. The human genome contains many pseudogenes similar to this locus.[provided by RefSeq, Feb 2011]