

## Product datasheet for **RC201976L4V**

### ACYP2 (NM\_138448) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ACYP2 (NM_138448) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ACYP2
Synonyms:	ACYM; ACYP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138448
ORF Size:	297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201976).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138448.2</a>
RefSeq Size:	1238 bp
RefSeq ORF:	300 bp



[View online »](#)

<b>Locus ID:</b>	98
<b>UniProt ID:</b>	<a href="#">P14621</a>
<b>Cytogenetics:</b>	2p16.2
<b>Protein Pathways:</b>	Pyruvate metabolism
<b>MW:</b>	11.1 kDa
<b>Gene Summary:</b>	<p>Acylphosphatase can hydrolyze the phosphoenzyme intermediate of different membrane pumps, particularly the Ca<sup>2+</sup>/Mg<sup>2+</sup>-ATPase from sarcoplasmic reticulum of skeletal muscle. Two isoenzymes have been isolated, called muscle acylphosphatase and erythrocyte acylphosphatase on the basis of their tissue localization. This gene encodes the muscle-type isoform (MT). An increase of the MT isoform is associated with muscle differentiation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016]</p>