

Product datasheet for **RC201986L4V**

CCN3 (NOV) (NM_002514) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CCN3 (NOV) (NM_002514) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCN3
Synonyms:	IBP-9; IGFBP-9; IGFBP9; NOV; NOVh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002514
ORF Size:	1071 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201986).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002514.2
RefSeq Size:	2601 bp
RefSeq ORF:	1074 bp
Locus ID:	4856
UniProt ID:	P48745
Cytogenetics:	8q24.12
Domains:	IB, tsp_1, VWC, CT, Cys_knot
Protein Families:	ES Cell Differentiation/IPS, Secreted Protein


[View online »](#)

MW: 39.2 kDa

Gene Summary: The protein encoded by this gene is a small secreted cysteine-rich protein and a member of the CCN family of regulatory proteins. CNN family proteins associate with the extracellular matrix and play an important role in cardiovascular and skeletal development, fibrosis and cancer development. [provided by RefSeq, Feb 2009]