


Product datasheet for RC201987L1

NEUROD1 (NM_002500) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NEUROD1 (NM_002500) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	NEUROD1
Synonyms:	BETA2; BHF-1; bHLHa3; MODY6; NEUROD; T2D
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201987).
Restriction Sites:	Sgfl-Mlul
Cloning Scheme:	
	Cloning sites used for ORF Shuttling:
	<i>Sgf I</i> ORF <i>Mlu I</i> GCG ATC GC <mark>ATG // NNŇ</mark> ACG CGT



ACCN: ORF Size: NM_002500 1068 bp

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


	NEURO	D1 (NM_002500) Human Tagged Lenti ORF Clone – RC201987L1
OTI Disclaimer:		Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <u>custsupport@origene.com</u> or by calling 301.340.3188 option 3 for pricing and delivery.
		The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:		This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:		The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution M	ethod:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:		<u>NM 002500.2</u>
RefSeq Size:		3002 bp
RefSeq ORF:		1071 bp
Locus ID:		4760
UniProt ID:		<u>Q13562</u>
Cytogenetics:		2q31.3
Domains:		HLH
Protein Families:		Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transcription Factors
Protein Pathways	:	Maturity onset diabetes of the young
MW:		39.9 kDa

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Gene Summary:This gene encodes a member of the NeuroD family of basic helix-loop-helix (bHLH)
transcription factors. The protein forms heterodimers with other bHLH proteins and activates
transcription of genes that contain a specific DNA sequence known as the E-box. It regulates
expression of the insulin gene, and mutations in this gene result in type II diabetes mellitus.
[provided by RefSeq, Jul 2008]

Product images:



Circular map for RC201987L1

Double digestion of RC201987L1 using Sgfl and Mlul

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US