

### Product datasheet for RC201989L2

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## EPCAM (NM\_002354) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: EPCAM (NM\_002354) Human Tagged Lenti ORF Clone

Tag: mGFP
Symbol: EPCAM

Synonyms: DIAR5; EGP-2; EGP40; EGP314; ESA; HNPCC8; KS1/4; KSA; M4S1; MIC18; MK-1; TACSTD1;

TROP1

Mammalian Cell

Selection:

None

Vector:pLenti-C-mGFP (PS100071)E. coli Selection:Chloramphenicol (34 ug/mL)

**ORF Nucleotide** 


The ORF insert of this clone is exactly the same as(RC201989).

Sequence:

**Restriction Sites:** Sgfl-Notl

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM 002354

ORF Size: 942 bp





#### EPCAM (NM\_002354) Human Tagged Lenti ORF Clone - RC201989L2

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 002354.1</u>, <u>NP 002345.1</u>

 RefSeq Size:
 1528 bp

 RefSeq ORF:
 945 bp

 Locus ID:
 4072

 UniProt ID:
 P16422

Cytogenetics: 2p21

Domains: thyroglobulin\_1

**Protein Families:** ES Cell Differentiation/IPS, Transmembrane

**MW:** 34.9 kDa


**Gene Summary:** This gene encodes a carcinoma-associated antigen and is a member of a family that includes

at least two type I membrane proteins. This antigen is expressed on most normal epithelial cells and gastrointestinal carcinomas and functions as a homotypic calcium-independent cell adhesion molecule. The antigen is being used as a target for immunotherapy treatment of human carcinomas. Mutations in this gene result in congenital tufting enteropathy. [provided

by RefSeq, Dec 2008]



# **Product images:**



Circular map for RC201989L2