

Product datasheet for RC201989L2V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

EPCAM (NM_002354) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: EPCAM (NM_002354) Human Tagged ORF Clone Lentiviral Particle

Symbol: EPCAM

Synonyms: DIAR5; EGP-2; EGP40; EGP314; ESA; HNPCC8; KS1/4; KSA; M4S1; MIC18; MK-1; TACSTD1;

TROP1

Mammalian Cell

Selection:

None

Vector: pLenti-C-mGFP (PS100071)

Tag: mGFP

ACCN: NM_002354

ORF Size: 942 bp

ORF Nucleotide

Sequence:

The ORF insert of this clone is exactly the same as(RC201989).

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 002354.1</u>, <u>NP 002345.1</u>

 RefSeq Size:
 1528 bp

 RefSeq ORF:
 945 bp

 Locus ID:
 4072

 UniProt ID:
 P16422

 Cytogenetics:
 2p21

Domains: thyroglobulin_1





EPCAM (NM_002354) Human Tagged ORF Clone Lentiviral Particle - RC201989L2V

Protein Families: ES Cell Differentiation/IPS, Transmembrane


MW: 34.9 kDa

Gene Summary: This gene encodes a carcinoma-associated antigen and is a member of a family that includes

at least two type I membrane proteins. This antigen is expressed on most normal epithelial cells and gastrointestinal carcinomas and functions as a homotypic calcium-independent cell adhesion molecule. The antigen is being used as a target for immunotherapy treatment of human carcinomas. Mutations in this gene result in congenital tufting enteropathy. [provided

by RefSeq, Dec 2008]

Product images:



[RC201989L2] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC201989L2V particle to overexpress human EPCAM-mGFP fusion protein.