

Product datasheet for RC202000L2V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD9 (NM_001769) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: CD9 (NM_001769) Human Tagged ORF Clone Lentiviral Particle

Symbol: CD9

Synonyms: BTCC-1; DRAP-27; MIC3; MRP-1; TSPAN-29; TSPAN29

Mammalian Cell

Selection:

None

Vector: pLenti-C-mGFP (PS100071)

Tag: mGFP

ACCN: NM_001769

ORF Size: 684 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC202000).

Sequence:

Cytogenetics:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 001769.2</u>

 RefSeq Size:
 1321 bp

 RefSeq ORF:
 687 bp

 Locus ID:
 928

 UniProt ID:
 P21926

Domains: transmembrane4

12p13.31



CD9 (NM_001769) Human Tagged ORF Clone Lentiviral Particle - RC202000L2V

Protein Families: Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS,

Transmembrane

Protein Pathways: Hematopoietic cell lineage

MW: 25.4 kDa

Gene Summary: This gene encodes a member of the transmembrane 4 superfamily, also known as the

tetraspanin family. Tetraspanins are cell surface glycoproteins with four transmembrane domains that form multimeric complexes with other cell surface proteins. The encoded protein functions in many cellular processes including differentiation, adhesion, and signal transduction, and expression of this gene plays a critical role in the suppression of cancer cell

motility and metastasis. [provided by RefSeq, Jan 2011]