

Product datasheet for **RC202005L4V**

C1orf156 (METTL18) (NM_033418) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C1orf156
Synonyms:	AsTP2; C1orf156; HPM1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_033418
ORF Size:	1116 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC202005).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_033418.1
RefSeq Size:	1487 bp
RefSeq ORF:	1119 bp
Locus ID:	92342
UniProt ID:	O95568
Cytogenetics:	1q24.2



MW: 42.1 kDa

Gene Summary: Probable histidine methyltransferase.[UniProtKB/Swiss-Prot Function]