

## Product datasheet for **RC202010L1V**

### Pea3 (ETV4) (NM\_001079675) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pea3 (ETV4) (NM_001079675) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Pea3
Synonyms:	E1A-F; E1AF; PEA3; PEAS3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001079675
ORF Size:	1452 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202010).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001079675.1</a>
RefSeq Size:	2469 bp
RefSeq ORF:	1455 bp
Locus ID:	2118
UniProt ID:	<a href="#">P43268</a>
Cytogenetics:	17q21.31
Protein Families:	Druggable Genome, Transcription Factors
MW:	53.9 kDa



[View online »](#)

**Gene Summary:** Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'. [UniProtKB/Swiss-Prot Function]