

Product datasheet for **RC202017L4V**

APPBP2 (NM_006380) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	APPBP2 (NM_006380) Human Tagged ORF Clone Lentiviral Particle
Symbol:	APPBP2
Synonyms:	APP-BP2; HS.84084; PAT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006380
ORF Size:	1755 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202017).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006380.2
RefSeq Size:	6468 bp
RefSeq ORF:	1758 bp
Locus ID:	10513
UniProt ID:	Q92624
Cytogenetics:	17q23.2
Domains:	TPR
Protein Families:	Druggable Genome



[View online »](#)

MW: 66.7 kDa

Gene Summary: The protein encoded by this gene interacts with microtubules and is functionally associated with beta-amyloid precursor protein transport and/or processing. The beta-amyloid precursor protein is a cell surface protein with signal-transducing properties, and it is thought to play a role in the pathogenesis of Alzheimer's disease. The encoded protein may be involved in regulating cell death. This gene has been found to be highly expressed in breast cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]