

Product datasheet for **RC202037L3V**

TPST1 (NM_003596) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TPST1 (NM_003596) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TPST1
Synonyms:	TANGO13A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003596
ORF Size:	1110 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202037).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003596.2
RefSeq Size:	2062 bp
RefSeq ORF:	1113 bp
Locus ID:	8460
UniProt ID:	O60507
Cytogenetics:	7q11.21
Domains:	Sulfotransfer
MW:	42.2 kDa


[View online »](#)

Gene Summary:

Catalyzes the O-sulfation of tyrosine residues within acidic motifs of polypeptides, using 3'-phosphoadenylyl sulfate (PAPS) as cosubstrate.[UniProtKB/Swiss-Prot Function]