

Product datasheet for **RC202048L4V**

TM9SF1 (NM_006405) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	TM9SF1
Synonyms:	HMP70; MP70
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006405
ORF Size:	1818 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC202048).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_006405.5
RefSeq Size:	2517 bp
RefSeq ORF:	1821 bp
Locus ID:	10548
UniProt ID:	O15321
Cytogenetics:	14q12



Domains:	EMP70
Protein Families:	Transmembrane
MW:	68.9 kDa
Gene Summary:	Plays an essential role in autophagy.[UniProtKB/Swiss-Prot Function]