

Product datasheet for **RC202050A1V**

Human RAIDD (CRADD) (NM_003805) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RAIDD (CRADD) (NM_003805) AAV Particle
Tag: Myc-DDK
Symbol: RAIDD
Synonyms: MRT34; RAIDD
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202050 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGCCAGAGACAAACAAGTACTCCGCTCACTTCGCCTGGAGCTGGGTGCAGAGGTATTGGTGGAGG
 GACTGGTTCTTCAGTACCTCTACCAGGAAGGAATCTTGACGGAACCATATTCAAGAAATCAATGCTCA
 AACCACAGGCCCTCCGAAAACAATGCTCCTGCTGGATATCCTACCTCCAGGGGCCCTAAAGCATTTGAT
 ACATTCCTAGATTCCCTACAGGAGTTTCCCTGGGTGAGGAGAAGCTGAAGAAGGCAAGGAAGAGGCCA
 TGACCGACCTGCCTGCAGGTGACAGATTGACTGGGATCCCTCGCACATCCTCAACAGCTCCCCATCAGA
 CCGGCAGATTAAACCAGCTGGCCAGAGGCTGGGCCCTGAGTGGGAGCCCATGGTGTCTCTGGGACTG
 TCCCAGACGGATATCTACCGCTGTAAGGCCAACCAACCCCAACAGTGCAGTGCAGGTGGTGGAGGCCT
 TCATCCGTTGGCGGCAGCGCTTCGGGAAGCAGGCCACCTTCAGAGCCTGCACAACGGGCTGCGGGCTGT
 GGAGGTGGACCCCTCGCTGCTCCTGCACATGTTGGAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202050 protein sequence
 Red=Cloning site Green=Tags(s)

MEARDKQVLRSLRLELGAEVLEGLVLQYLYQEGILTENHIQEINAQTTGLRKTMLLLDILPSRGPKA
 FDTFLDSLQEFPPWREKLKAREEAMTDLPAAGDRLTGIPSHILNSSPSDRQINQLAQRLGPEWEP
 MVLSLGLSQTDIYRCKANHPHNVSQVVEAFIRWRQRFQKQATFQSLHNGLRAVEVDP
 SLLHMLE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_003805
ORF Size:	597 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003805.3
RefSeq Size:	1201 bp
RefSeq ORF:	600 bp
Locus ID:	8738
UniProt ID:	P78560
Cytogenetics:	12q22
MW:	22.7 kDa