

Product datasheet for **RC202051A1V**

Human D Box Binding Protein (DBP) (NM_001352) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human D Box Binding Protein (DBP) (NM_001352) AAV Particle
Tag:	Myc-DDK
Symbol:	D Box Binding Protein
Synonyms:	DABP; taxREB302
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202051 representing NM_001352 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGCGCCTGTGAGCGACAGGACCCCGGCCCTCTGCTGCTGGGCGGCCCGGCCGGACACCCCTG
GCGGGGAGCGCTGCTTGGGTTGCGGAGCCTTCTGCAGGGGACCAGCAAGCCAAAGAGCCGGCCAGCTG
TCTCTGAAGGAAAAGGAGCGCAAGGCGGCCCTGCCTGCAGCCACAACCCCTGGGCCAGGCCTGGAGACT
GCGGGCCCGCGGATGCCCGGCTGGGCGAGTGGTGGGCGGAGGGTCCCCGCGGGGGCGCCCGGGCCGG
TGCCCGCCCGGGTCTGTTGGCGCCACTGCTGTGGGAGCGCACGCTGCCGTTTCGGCGATGTGGAGTACGT
AGACCTGGACGCCTTCTGCTGGAGCACGGGCTCCCGCCAGCCCGCCCGCCCGGTCGGCCGTCGCCG
GAGCCGTCGCCCGCGCGGACGCCCGCACCCCTCCCGAGGCGGGTTCGTGCGGCTCGGCTTCCCCCGCT
CCTCTCTGGGACGCCCCCGCCGGGCTGCCCTCGGGACCGCCAGCGGCCACCGCGCAGGCCTGACCTC
TCGGGACACACCCAGCCCTGTGGACCCAGACACCGTGGAGGTGTTGATGACCTTTGAACCCGACCCAGCT
GATCTTGCCCTATCAAGCATTCTGGCCACGAGACCTTTGACCCTCGAAGACATCGCTTCTCAGAAGAGG
AACTTAAGCCCCAGCCAATCATGAAGAAGGCAAGAAAAATCCAGGTGCCGGAGGAGCAGAAGGATGAGAA
ATACTGGAGCCGGCGGTACAAGAACAACGAGGCAGCCAAAGCGGTCCCGTGACGCCCGCGGCTCAAGGAG
AACCAGATATCGGTGCGGGCGGCCTTCTGGAGAAGGAGAACGCCCTGCTGCGGCAGGAAGTTGTGGCCG
TGCGCCAGGAGCTGTCCCACTACCGCGCGTGTGTCCCGATACCAAGGCCAGCACGGGGCCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202051 representing NM_001352 Red=Cloning site Green=Tags(s) MARPVSDRTPAPLLLGGPAGTPPGGGALLGLRSLLQGTSKPKEPASCLLKEKERKAALPAATTPGPGLET AGPADAPAGAVVGGGSPRGRPGVPAPGLLAPLLWERTLPFGDVEYVDLDAFLLEHGLPPSPPPGGPSP EPSPARTPAPSPGPGSCGSASPRSSPGHAPARAALGTASGHRAGLTSRDTPSPVDPDTVEVLMTFEPDPA DLALSSIPGHETFDPRRHRFSEEELKPQPIIMKKARKIQVPEEQKDEKYWSRRYKNNEAAKRSRDARRLKE NQISVRAAFLEKENALLRQEVVAVRQELSHYRAVL SRYQAQHGAL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001352
ORF Size:	975 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001352.2</u>
RefSeq Size:	1478 bp
RefSeq ORF:	978 bp
Locus ID:	1628
UniProt ID:	<u>Q10586</u>
Cytogenetics:	19q13.33
MW:	34.2 kDa