

Product datasheet for **RC202063L4V**

G0S2 (NM_015714) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	G0S2 (NM_015714) Human Tagged ORF Clone Lentiviral Particle
Symbol:	G0S2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015714
ORF Size:	309 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202063).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015714.2
RefSeq Size:	978 bp
RefSeq ORF:	312 bp
Locus ID:	50486
UniProt ID:	P27469
Cytogenetics:	1q32.2
Protein Families:	Transmembrane
MW:	11.3 kDa



[View online »](#)

Gene Summary:

Promotes apoptosis by binding to BCL2, hence preventing the formation of protective BCL2-BAX heterodimers.[UniProtKB/Swiss-Prot Function]