

Product datasheet for RC202064L2

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SIAT4A (ST3GAL1) (NM_173344) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: SIAT4A (ST3GAL1) (NM 173344) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: SIAT4A

Synonyms: 1; Gal-NAc6S; SIAT4A; SIATFL; ST3GalA; ST3GalA.1; ST3GalIA; ST3O

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC202064).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_173344

ORF Size: 1020 bp



SIAT4A (ST3GAL1) (NM_173344) Human Tagged Lenti ORF Clone - RC202064L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 173344.1</u>

 RefSeq Size:
 6768 bp

 RefSeq ORF:
 1023 bp

 Locus ID:
 6482

 UniProt ID:
 011201

Cytogenetics: 8q24.22

Protein Families: Secreted Protein, Transmembrane

Protein Pathways: Glycosphingolipid biosynthesis - ganglio series, Glycosphingolipid biosynthesis - globo series,

Keratan sulfate biosynthesis, Metabolic pathways, O-Glycan biosynthesis

MW: 39.1 kDa


Gene Summary: The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of

sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein is normally found in the Golgi but can be proteolytically processed to a soluble form. Correct glycosylation of the encoded protein may be critical to its sialyltransferase activity. This protein, which is a member of glycosyltransferase family 29, can use the same acceptor substrates as does sialyltransferase 4B. Two transcript variants encoding the same protein have been found for this gene. Other transcript variants may exist, but have not been fully

characterized yet. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC202064L2