

Product datasheet for **RC202071A1V**

Human Noxa (PMAIP1) (NM_021127) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Noxa (PMAIP1) (NM_021127) AAV Particle
Tag: Myc-DDK
Symbol: Noxa
Synonyms: APR; NOXA
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202071 representing NM_021127
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCTGGGAAGAAGGCGCGCAAGAACGCTCAACCGAGCCCCGCGGGCTCCAGCAGAGCTGGAAGTCG
AGTGTGCTACTCAACTCAGGAGATTTGGAGACAACTGAACTTCCGGCAGAACTTCTGAATCTGATATC
CAAACCTTCTGCTCAGGAACC

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202071 representing NM_021127
 Red=Cloning site Green=Tags(s)
 MPGKKARKNAQPSPARAPAELEVECATQLRRFGDKLNFRQKLLNLISKLFCSGT
TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_021127
ORF Size: 162 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_021127.1](#)

RefSeq Size: 1885 bp

RefSeq ORF: 165 bp

Locus ID: 5366

UniProt ID: [Q13794](#)

Cytogenetics: 18q21.32

MW: 5.8 kDa