

Product datasheet for **RC202076L4V**

MLC1SA (MYL6B) (NM_002475) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MLC1SA (MYL6B) (NM_002475) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MLC1SA
Synonyms:	MLC1SA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002475
ORF Size:	624 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202076).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002475.2
RefSeq Size:	1008 bp
RefSeq ORF:	627 bp
Locus ID:	140465
UniProt ID:	P14649
Cytogenetics:	12q13.2
Domains:	EFh
Protein Pathways:	Vascular smooth muscle contraction



[View online »](#)

MW: 22.8 kDa

Gene Summary: Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain expressed in both slow-twitch skeletal muscle and in nonmuscle tissue. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010]