

## Product datasheet for **RC202085L1V**

### **PAI1 (SERPINE1) (NM\_000602) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PAI1 (SERPINE1) (NM_000602) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PAI1
Synonyms:	PAI; PAI-1; PAI1; PLANH1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000602
ORF Size:	1206 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202085).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000602.1</a>
RefSeq Size:	3207 bp
RefSeq ORF:	1209 bp
Locus ID:	5054
UniProt ID:	<a href="#">P05121</a>
Cytogenetics:	7q22.1
Domains:	SERPIN
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**Protein Pathways:** Complement and coagulation cascades, p53 signaling pathway

**MW:** 45.1 kDa

**Gene Summary:** This gene encodes a member of the serine proteinase inhibitor (serpin) superfamily. This member is the principal inhibitor of tissue plasminogen activator (tPA) and urokinase (uPA), and hence is an inhibitor of fibrinolysis. The protein also functions as a component of innate antiviral immunity. Defects in this gene are the cause of plasminogen activator inhibitor-1 deficiency (PAI-1 deficiency), and high concentrations of the gene product are associated with thrombophilia. [provided by RefSeq, Aug 2020]