

Product datasheet for **RC202090A1V**

Human RBBP9 (NM_006606) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Human RBBP9 (NM_006606) AAV Particle
 Tag: Myc-DDK
 Symbol: RBBP9
 Synonyms: BOG; RBBP10
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >RC202090 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTTCTCCTAGCAAGGCAGTGATTGTTCCCGGGAACGGAGGCGGGGATGTGACCACCCACGGCTGGT
 ATGGCTGGGTGAAAAAGGAGCTGGAGAAGATACCTGGTTTCCAGTGTTTGGCTAAAAACATGCCCGACCC
 AATTACAGCACGAGAGAGCATCTGGCTGCCCTTCATGGAGACAGAGCTGCACTGTGATGAGAAGACTATC
 ATCATTGGCCACAGTTCTGGGGCCATCGCGCCATGAGGTATGCAGAAACACATCGAGTATATGCTATTG
 TATTAGTGTCTGCGTACACATCAGACTTGGGGGATGAAATGAGCGTGCAAGTGGATACTTCACCCGCC
 CTGGCAGTGGGAGAAGATCAAGGCCAACTGCCCTTACATTGTGCAGTTTGGCTCTACTGACGACCCGTTT
 CTTCCCTGGAAGGAACAACAAGAAGTGCCGATAGGTTGGAACCAAATTGCACAAATCACTGACTGTG
 GCCACTTTCAGAACACAGAGTTTCATGAAGTACTGTTGTAAGTCTTGTCTGAAAGTACCAGCA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202090 protein sequence
 Red=Cloning site Green=Tags(s)

MASPSKAVIVPGNGGGDVTTHGWYGWVKELEKIPGFQCLAKNMPDPITARESIWLPFMEtelHCDEKTI
 IIGHSSGAIAAMRYAETHRVYAIVLVSAYTSDLGDNERASGYFTRPWQWEKIKANCPYIVQFGSTDDPF
 LPWKEQQEVADRLETKLHKFTDCGHFQNTFHELITVVKSLKVPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
 Serotype: AAV-2



[View online »](#)

ACCN:	NM_006606
ORF Size:	558 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006606.2
RefSeq Size:	3871 bp
RefSeq ORF:	561 bp
Locus ID:	10741
UniProt ID:	O75884
Cytogenetics:	20p11.23
MW:	21 kDa