

Product datasheet for **RC202098A1V**

Human TFPI (NM_001032281) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TFPI (NM_001032281) AAV Particle
Tag:	Myc-DDK
Symbol:	TFPI
Synonyms:	EPI; LACI; TFI; TFPI1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202098 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGATTTACACAATGAAGAAAGTACATGCACTTTGGGCTTCTGTATGCCTGCTGCTTAATCTTGCCCTG
CCCCTCTTAATGCTGATTCTGAGGAAGATGAAGAACACACAATTATCACAGATACGGAGTTGCCACCACT
GAACTTATGCATTCATTTTGTGCATTCAAGGCGGATGATGGCCCATGTAAGCAATCATGAAAAGATTT
TTCTTCAATATTTTCACTCGACAGTGCGAAGAATTTATATATGGGGGATGTGAAGGAAATCAGAATCGAT
TTGAAAGTCTGGAAGAGTGCAAAAAAATGTGTACAAGAGATAATGCAACAGGATTATAAGACAACATT
GCAACAAGAAAAGCCAGATTTCTGCTTTTTGGAAGAAGATCCTGGAATATGTCGAGGTTATATTACCAGG
TATTTTATAACAATCAGACAAAACAGTGTGAACGTTTCAAGTATGGTGGATGCCTGGGCAATATGAACA
ATTTTGAGACACTGGAAGAATGCAAGAACATTTGTGAAGATGGTCCGAATGGTTTCCAGGTGGATAATTA
TGGAACCCAGCTCAATGCTGTGAATAACTCCCTGACTCCGCAATCAACCAAGTTCCAGCCTTTTTGTT
ACAAAAGAAGGAACAAATGATGGTTGGAAGAATGCGGCTCATATTTACCAAGTCTTTCTGAACGCCTTCT
GCATTCATGCATCCATGTTCTTTCTAGGATTGGATAGCATTTTCATGCCTATGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202098 protein sequence Red=Cloning site Green=Tags(s) MIYTMKKVHALWASVCLLLNLAPAPLNADSEEDDEHTIITDTELPPLKLMHSFCAFKADDGPCKAIMKRF FFNIFTRQCEEFYGGCEGNQNRFSLEECKMCTRDNANRIIKTTLQKEKPDFCFLEEDPGICRGYITR YFYNNQTKQCERFKYGGCLGNMNNFETLEECKNICEDGPNQFQVDNYGTQLNAVNNSLTPQSTKVPRLFV TKEGTNDGWKNAAHYQVFLNAFCIHASMFFLGLDNISCLC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001032281
ORF Size:	753 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001032281.2, NP_001027452.1</u>
RefSeq Size:	1166 bp
RefSeq ORF:	756 bp
Locus ID:	7035
UniProt ID:	<u>P10646</u>
Cytogenetics:	2q32.1
MW:	28.7 kDa