

Product datasheet for **RC202102A1V**

Human IRF2 (NM_002199) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human IRF2 (NM_002199) AAV Particle
Tag:	Myc-DDK
Symbol:	IRF2
Synonyms:	IRF-2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202102 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCCGGTGGAAGGATGCGCATGCGCCGTGGCTGGAGGAGCAGATAAACTCCAACACGATCCCGGGGC
TCAAGTGGCTTAACAAGGAAAAGAAGATTTTCAGATCCCCTGGATGCATGCGGCTAGACATGGGTGGGA
TGTGGAAAAAGATGCACCACTCTTTAGAACTGGGCAATCCATACAGGAAAGCATCAACCAGGAGTAGAT
AAACCTGATCCCAAAACATGGAAGGCGAATTCAGATGCGCCATGAATTCCTTGCCTGATATTGAAGAAG
TCAAGGATAAAAGCATAAAGAAAGGAAATAATGCCTTCAGGGTCTACCGAATGCTGCCCTATCAGAACG
GCCTTCTAAGAAAGGAAAGAAACCAAGACAGAAAAAGACAAAGTTAAGCACATCAAGCAAGAACCA
GTTGAGTCATCTCTGGGGCTTAGTAATGGAGTAAGTGATCTTTCTCCTGAGTATGCGGTCTGACTTCAA
CTATAAAAAATGAAGTGGATAGTACGGTGAACATCATAGTTGTAGGACAGTCCCATCTGGACAGCAACAT
TGAGAATCAAGAGATTGTCACCAATCCGCCAGACATTTGCCAAGTTGTAGAGGTGACCACTGAGAGCGAC
GAGCAGCCGTCAGCATGAGCGAGCTCTACCCTCTGCAGATCTCCCCGTGTCTTCTATGCAGAAAGCG
AAACGACTGATAGTGTGCCAGCGATGAAGAGAGTGCCGAGGGGCGGCCACACTGGCGGAAGAGGAATAT
TGAAGGCAACAGTACCTCAGCAACATGGGACTCGAGGCTCTACCTGCTGCCGGCATGGCGTCTTTC
GTCACCTTCCAACAAACCGGACCTCCAGGTCACCATCAAAGAGGAGAGCAATCCGGTGCCTTACAACAGCT
CCTGGCCCCCTTTTCAAGACCTCCCCCTTTCTTCTCCATGACCCAGCATCCAGCAGCAGTCCGGCCAGA
CCGGGAGACCCGGGCCAGCGTCATCAAGAAAACATCGGATATCACCCAGGCCCGCTCAAGAGCTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202102 protein sequence Red=Cloning site Green=Tags(s) MPVERMRMPWLEEQINSNTIPGLKWLNKEKKIFQIPWMHAARHGWDVEKDAPLFRNWAHTGKHQPGVD KPDPKTWKANFRCAMNSLPDIEEVKDKSIKKGNNAFRVYRMLPLSERPSKKGKKPKTEKEDKVKKHIKQEP VESSLGLSNGVSDLSPEYAVLTSTIKNEVDSTVNIIVVGQSHLDSNIENQEIVTNPPDICQVVEVTESD EQPVSMSELYPLQISPVSSYAESETTDSVPSDEESAEGRPHWRKRNIIEGKQYLSNMGTRGSYLLPGMASF VTSNKPDLQVTIKEESNPVPYNSSWPPFQDLPLSSMTPASSSSRPDRETRASVIKKTSDITQARVKSC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002199
ORF Size:	1047 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_002199.2
RefSeq Size:	2302 bp
RefSeq ORF:	1050 bp
Locus ID:	3660
UniProt ID:	P14316
Cytogenetics:	4q35.1
MW:	39.4 kDa