

Product datasheet for **RC202108L1V**

Selenium Binding Protein 1 (SELENBP1) (NM_003944) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Selenium Binding Protein 1 (SELENBP1) (NM_003944) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Selenium Binding Protein 1
Synonyms:	EHMTO; HEL-S-134P; hSBP; LPSB; MTO; SBP56; SP56
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003944
ORF Size:	1416 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202108).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003944.2
RefSeq Size:	1768 bp
RefSeq ORF:	1419 bp
Locus ID:	8991
UniProt ID:	Q13228
Cytogenetics:	1q21.3
MW:	52.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the selenium-binding protein family. Selenium is an essential nutrient that exhibits potent anticarcinogenic properties, and deficiency of selenium may cause certain neurologic diseases. The effects of selenium in preventing cancer and neurologic diseases may be mediated by selenium-binding proteins, and decreased expression of this gene may be associated with several types of cancer. The encoded protein may play a selenium-dependent role in ubiquitination/deubiquitination-mediated protein degradation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Apr 2012]