

Product datasheet for **RC202109A1V**

Human PTTG1IP (NM_004339) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PTTG1IP (NM_004339) AAV Particle
Tag: Myc-DDK
Symbol: PTTG1IP
Synonyms: C21orf1; C21orf3; PBF
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202109 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGCCCGGAGTGGCCCGGGCCGACGCCGTAAGGTTGCGCCTCGGTGGCGCCGCGTCTCTCC
TGCTGCTCATCCCGGTGGCCCGCGCAGGAGCCTCCCGGAGCTGCTTGTCTCAGAACAACAAAAAC
CTGTGAAGAGTGCCTGAAGAAGTCTCTGTCTTTGGTGCAACACTAACAGGCTTGTCTGGACTACCA
GTTACAAGCGCTTGGCACCAGGCTTCCCTTTGTAATTGAGCTCTGCACGCTGGGAGTTTGTGGGTGA
ACTTTGAGGCGCTGATCATCACCATGTCGGTAGTCGGGGGAACCTCCTCTGGGCATTGCCATCTGCTG
CTGTGCTGCTGCAGGAGGAAGAGGAGCGGAAGCCGGACAGGAGTGAGGAGAAGGCCATGCGTGAGCGG
GAGGAGAGGCGGATACGGCAGGAGGAACGGAGAGCAGAGATGAAGACAAGACATGATGAAATCAGAAAAA
AATATGGCCTGTTTAAAGAAGAAAACCGTATGCTAGATTGAAAAAAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202109 protein sequence
Red=Cloning site Green=Tags(s)

MAPGVARGPTYWRLRLGGAALLLLLIPVAAAEPPGAACSQNTNKTCEECLKNVSLWCNTNKACLDYP
VTSVLPPASLCKLSSARWGVWNFEALITMSVVGTTLLGLAICCCCCRRKRSRKPDRSEEKAMRER
EERRIRQEERRAEMKTRHDEIRKKYGLFKEENPYARFENN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



View online »

ACCN:	NM_004339
ORF Size:	540 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_004339.2</u>
RefSeq Size:	2843 bp
RefSeq ORF:	543 bp
Locus ID:	754
UniProt ID:	<u>P53801</u>
Cytogenetics:	21q22.3
MW:	20.3 kDa