

## Product datasheet for **RC202119L2V**

### **PDCD5 (NM\_004708) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PDCD5 (NM_004708) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PDCD5
Synonyms:	TFAR19
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004708
ORF Size:	375 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202119).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004708.2</a>
RefSeq Size:	604 bp
RefSeq ORF:	378 bp
Locus ID:	9141
UniProt ID:	<a href="#">O14737</a>
Cytogenetics:	19q13.11
Domains:	dsDNA_bind
MW:	14.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that is upregulated during apoptosis where it translocates rapidly from the cytoplasm to the nucleus. The encoded protein may be an important regulator of K(lysine) acetyltransferase 5 (a protein involved in transcription, DNA damage response and cell cycle control) by inhibiting its proteasome-dependent degradation. Pseudogenes have been identified on chromosomes 5 and 12 [provided by RefSeq, Dec 2010]