

## Product datasheet for **RC202125L3V**

### TCIRG1 (NM\_006019) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TCIRG1 (NM_006019) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TCIRG1
Synonyms:	a3; Atp6i; ATP6N1C; ATP6V0A3; OC-116kDa; OC116; OPTB1; Stv1; TIRC7; Vph1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006019
ORF Size:	2490 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202125).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006019.2</a>
RefSeq Size:	2727 bp
RefSeq ORF:	2493 bp
Locus ID:	10312
UniProt ID:	<a href="#">Q13488</a>
Cytogenetics:	11q13.2
Domains:	V_ATPase_sub_a
Protein Families:	Transmembrane



[View online »](#)

<b>Protein Pathways:</b>	Epithelial cell signaling in Helicobacter pylori infection, Lysosome, Metabolic pathways, Oxidative phosphorylation, Vibrio cholerae infection
<b>MW:</b>	93 kDa
<b>Gene Summary:</b>	This gene encodes a subunit of a large protein complex known as a vacuolar H <sup>+</sup> -ATPase (V-ATPase). The protein complex acts as a pump to move protons across the membrane. This movement of protons helps regulate the pH of cells and their surrounding environment. V-ATPase dependent organelle acidification is necessary for such intracellular processes as protein sorting, zymogen activation, and receptor-mediated endocytosis. V-ATPase is comprised of a cytosolic V1 domain and a transmembrane V0 domain. Alternative splicing results in multiple transcript variants. Mutations in this gene are associated with infantile malignant osteopetrosis. [provided by RefSeq, May 2017]