

## Product datasheet for **RC202126L1V**

### UCH37 (UHL5) (NM\_015984) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UCH37 (UHL5) (NM_015984) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UCH37
Synonyms:	CGI-70; INO80R; UCH-L5; UCH37
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_015984
ORF Size:	987 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202126).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015984.1</a>
RefSeq Size:	1728 bp
RefSeq ORF:	990 bp
Locus ID:	51377
UniProt ID:	<a href="#">Q9Y5K5</a>
Cytogenetics:	1q31.2
Domains:	Peptidase_C12
Protein Families:	Druggable Genome, Protease



[View online »](#)

**MW:** 37.4 kDa

**Gene Summary:** Protease that specifically cleaves 'Lys-48'-linked polyubiquitin chains. Deubiquitinating enzyme associated with the 19S regulatory subunit of the 26S proteasome. Putative regulatory component of the INO80 complex; however is inactive in the INO80 complex and is activated by a transient interaction of the INO80 complex with the proteasome via ADRM1. [UniProtKB/Swiss-Prot Function]