

## **Product datasheet for RC202135L1**

#### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

### ASAH1 (NM\_177924) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: ASAH1 (NM 177924) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: ASAH1

Synonyms: AC; ACDase; ASAH; PHP; PHP32; SMAPME

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC202135).

Sequence:

**Restriction Sites:** Ascl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_177924

ORF Size: 1185 bp





### ASAH1 (NM\_177924) Human Tagged Lenti ORF Clone - RC202135L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 177924.2</u>, <u>NP 808592.1</u>

RefSeq Size:2609 bpRefSeq ORF:1188 bp

Locus ID: 427

UniProt ID: Q13510

Cytogenetics: 8p22

**Protein Families:** Druggable Genome

**Protein Pathways:** Lysosome, Metabolic pathways, Sphingolipid metabolism

MW: 42.6 kDa


**Gene Summary:** This gene encodes a member of the acid ceramidase family of proteins. Alternative splicing

results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed. Processing of this preproprotein generates alpha and beta subunits that heterodimerize to form the mature lysosomal enzyme, which catalyzes the degradation of ceramide into sphingosine and free fatty acid. This enzyme is overexpressed in multiple human cancers and may play a role in cancer progression. Mutations in this gene are associated with the lysosomal storage disorder, Farber lipogranulomatosis, and a neuromuscular disorder, spinal muscular atrophy with progressive myoclonic epilepsy.


[provided by RefSeq, Oct 2015]



# **Product images:**



Circular map for RC202135L1



Double digestion of RC202135L1 using Sgfl and Mlul  $\,$